

Drug Cupboards

Controlled Drugs Cabinet

- Designed and manufactured in mild steel to meet "The Misuse of Drugs (Safe Custody) Regulations 1973" (except CD040, CD045, CD050 and CD055)
- Recommended by "The Home Office Drugs Inspectorate" for storage of Controlled Drugs
- Security can be enhanced by concealing inside a PC cabinet
- Single-point locking
- Fitted with a pick and drill resistant key conforming to BS3621
- Compatible with electronic key and keyless locking systems
- Can be fitted with a micro-switch, to allow connection to a volt-free contact, i.e. nurse call system/remote lighting
- Finished in a durable lead free epoxy/polyester powder coated paint
- Fixing holes provided in back
- Adjustable and removable shelves, (refer to table models)
- Modular design allows for easy interchange with wooden fitted storage systems, (HTM 63)
- w 210 x d 270 x h 300mm (other sizes are available)

CODE	DESCRIPTION	QTY.	PRICE
CD005	Controlled Drug Cupboard	1	£173.30

Drug Cupboard with MDS System

- Complies with BS2881:1989, (Security Level 1)
- Adjustable and removable shelving
- Unique 12v DC lighting/warning system which is CE marked
- Nurse call system ready (volt free contact)
- Five pin key retaining lock
- Modular design allows for easy interchange with wooden fitted storage systems, (HTM 63)
- Transformer casing adds 105mm to the height
- w 500 x d 300 x h 900mm (other sizes are available)

CODE	DESCRIPTION	QTY.	PRICE
PC075	Cupboard 4 Frame Single Door	1	£340.70



Drug Trolley

Monitored Dosage System Trolley

- A new range of lockable trolleys designed for the storage and dispensing of the popular Blister Packed Monitored Dosage System, (MDS)
- Push handles fitted to both ends incorporating integral buffers to protect the fabric of the building as well as assist manoeuvring the trolley
- Frames supported on either a fixed panel mounted to the inside rear of the trolley, or on adjustable slide out shelves
- Models available for the storage of 6 and 9 frames
- Fold over flap fitted to the door complete with integral bracket to support MDS frames whilst dispensing
- Removable and adjustable door shelves to accommodate bottles/packs.
- Plastic dispensing tray
- Plastic debris bin with bracket fitted to one side of the trolley
- Fitted with 125mm swivel castors
- Drawer available for the storage of loose items
- Can be secured to the fabric of the building using a security clamp MT/1/6F (see below)
- w 670 x d 530 x h 1095mm (other sizes are available)




CODE	DESCRIPTION	QTY.	PRICE
UD210	Monitored Dosage System Trolley 6 Frame	1	£640.00

Security Clamp

Security clamp for use with Monitored Dosage System Trolley.

- Secures drug trolley to the fabric of the building when not in use
- No additional keys required, as it utilises the trolley's own locking device
- The clamp fits in to a slot in the top of the door and is trapped when the trolley is locked
- Recommended fitting height approx. 100cm



CODE	DESCRIPTION	QTY.	PRICE
MT-1-6F	Security Clamp	1	£24.00

Full Brochure Available

Please telephone us on 01773 713713 or email us at sales@hcsuk.co.uk to request a Bristol Maid full range brochure.



Please allow 6 - 8 weeks for delivery on all Bristol Maid drug orders